



AutoMolec FluA/FluB/RSV
AutoMolec HPIV
AutoMolec HMPV
AutoMolec M. Pneumoniae/Adenovirus
AutoMolec SA/MRSA/KPN
AutoMolec SP/LP/MC
AutoMolec PAE/HIN/E.coli
AutoMolec M. Pneumoniae & Macrolides -Resistance
AutoMolec M. Pneumoniae/C. Pneumoniae/Adenovirus*
AutoMolec Victoria/Yamagata*
AutoMolec AIV H5/H7*
AutoMolec Flu A H1/H3*

AutoMolec Norovirus*

AutoMolec Rotavirus group A*

AutoMolec Salmonella/Shigella*

AutoMolec HPV 16/18
AutoMolec GBS
AutoMolec HSV 1/2
AutoMolec CT/NG/UU*
AutoMolec HPV 6/11*

AutoMolec CYP2C19

AutoMolec ALDH2

AutoMolec MTHFR/MTRR

AutoMolec HBV DNA Quant*

AutoMolec HCV RNA Quant*

AutoMolec Monkeypox *

*RUO

INNOVATION · MANUFACTURE · INTEGRATION · SERVICE

Autobio Diagnostics Co., Ltd.

offers more than 600 in vitro diagnostic products including CLIA(microplate based CLIA and magnetic particle based CLIA),ELISA,POCT(Point of Care Test),Microbiology and Biochemistry. As an ISO9001 and EN ISO13485 manufacturer, Autobio supplies high quality products through its well established sales network and is renowned as a reliable partner. For details please visit en.autobio.com.cn. Autobio Diagnostics Co., Ltd. | NO.87 Jingbei Yi Rd | National Eco & Tech Zone | Zheng zhou City | China 450016.

You can find us on:

Facebook: Autbio Diagnostics
 LinkedIn: Autbio Diagnostics
 YouTube Channel: Autbio Diagnostics
 Email: info@autbio-diagnostics.com
 WhatsApp: +86 185 9576 7835
 Tel: +86 371 6200 7036



Automatic molecular diagnostic system20240627e

Automated Nucleic Acid Purification and Real Time PCR System



en.autobio.com.cn




AutoMoloc 3000—Fully integrated, Completely automated, Entirely flexible

The AutoMolec 3000 is a fully integrated and automated molecular diagnostics analyzer, which integrates nucleic acid extraction and real-time PCR in a single workflow. 41 independent temperature controlling modules ensure true random access. The analyzer system makes your nucleic acid assessment more accurate, flexible and convenient.


2. Sample inlet

- ※ Random access
- ※ STAT inlet



3. Extraction

- ※ 20 mins for first sample
- ※ 1 more sample every 2 mins




1. Anti-contamination

- ※ Independant dual air control
- ※ Intelligent UV light
- ※ Triple anti-droplet system in pipetting
- ※ Sealing system in PCR
- ※ Special designed to prevent cross-contamination dropping





4. Reagents and consumables

- ※ 20 slots with 2-8℃ cooling
- ※ Supporting 20 test items on-board
- ※ On-the-fly change for reagents and consumables




5. PCR system

- ※ 41 independant control channels
- ※ Supporting both RNA and DNA
- ※ RNA 80 min, DNA 50 mins



Working Flow



Specification

Model	AutoMolec3000
Sample Loading Capacity	60 samples(circular rotor), priority lane for emergency samples.
Projects Analyzed Simultaneously	20 items at one time
Kit Slot	20, online order and replacement.
Throughput	200tests/8h
Test Mode	It supports independent temperature control, independent testing, on-call test and sample purification and amplification integrated.
Reagent Temperature	2℃~8℃
Weight	670 kg
Size	2400 mm×900 mm×1900 mm (L×W×H)