

Infectious Disease Control III– Level 2

Assignment of values (**LOT** 15C0701)

Instrument	Analyte	Unit	Infectious Disease Control III-Level 2	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus	HBsAg	IU/mL	0.46	0.32-0.59
	HBeAg	IU/mL	0.33	0.23-0.43
	Anti-HIV	S/CO	3.14	2.20-4.08
	HIV Ag/Ab combo	S/CO	2.89	2.02-3.76
	Anti-TP	S/CO	3.05	2.13-3.96
	Anti-HCV	S/CO	3.03	2.12-3.93

Instrument	Analyte	Unit	Infectious Disease Control III-Level 2	
			Target	Range
Auto Lumo A1000	HBsAg	IU/mL	0.41	0.29-0.54
	HBeAg	IU/mL	0.30	0.21-0.38
	Anti-HIV	S/CO	2.83	1.98-3.67
	HIV Ag/Ab combo	S/CO	2.91	2.04-3.78
	Anti-TP	S/CO	2.73	1.91-3.55
	Anti-HCV	S/CO	2.84	1.99-3.70

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.