

## Infectious Disease Control III– Level 3

Assignment of values ( **LOT** 15C0701 )

Instrument	Analyte	Unit	Infectious Disease Control III-Level 3	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus	HBsAg	IU/mL	3.72	2.60-4.84
	HBeAg	IU/mL	3.75	2.62-4.87
	Anti-HIV	S/CO	6.21	4.34-8.07
	HIV Ag/Ab combo	S/CO	6.16	4.31-8.00
	Anti-TP	S/CO	5.93	4.15-7.71
	Anti-HCV	S/CO	5.76	4.03-7.48

Instrument	Analyte	Unit	Infectious Disease Control III-Level 3	
			Target	Range
Auto Lumo A1000	HBsAg	IU/mL	3.97	2.78-5.16
	HBeAg	IU/mL	3.58	2.51-4.65
	Anti-HIV	S/CO	5.58	3.91-7.26
	HIV Ag/Ab combo	S/CO	5.87	4.11-7.63
	Anti-TP	S/CO	5.32	3.72-6.92
	Anti-HCV	S/CO	5.39	3.77-7.01

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.