

Infectious Disease Control IV – Level 2

Assignment of values (**LOT** 15C0701)

Instrument	Analyte	Unit	Infectious Disease Control IV-Level 2	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus	Anti-HBs	mIU/mL	31.25	21.87-40.62
	Anti-HBe	PEI U/mL	0.50	0.35-0.65
	Anti-HBc	PEI U/mL	2.29	1.61-2.98
	HIV Ag/Ab combo	S/CO	3.01	2.11-3.92

Instrument	Analyte	Unit	Infectious Disease Control IV-Level 2	
			Target	Range
Auto Lumo A1000	Anti-HBs	mIU/mL	28.75	20.12-37.37
	Anti-HBe	PEI U/mL	0.53	0.37-0.69
	Anti-HBc	PEI U/mL	2.41	1.69-3.13
	HIV Ag/Ab combo	S/CO	2.83	1.98-3.68

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.