

Infectious Disease Control III– Level 4

Assignment of values(LOT 22K0701)

Instrument	Analyte	Unit	Infectious Disease Control III-Level 4	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus	HBsAg	IU/mL	37.30	26.11-48.48
	HBeAg	IU/mL	29.57	20.70-38.44
	Anti-TP	S/CO	12.31	8.62-16.01
	Anti-HCV	S/CO	13.11	9.18-17.04
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	12.19	8.53-15.84
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	12.45	8.71-16.18
Auto Lumo A1000	HBsAg	IU/mL	44.23	30.96-57.49
	HBeAg	IU/mL	30.56	21.39-39.72
	Anti-TP	S/CO	11.56	8.09-15.02
	Anti-HCV	S/CO	13.21	9.25-17.18
	Anti-HIV (Anti-HIV CLIA Microparticles)	S/CO	11.34	7.94-14.74
	Anti-HIV (HIV Ag/Ab Combo CLIA Microparticles)	S/CO	11.77	8.24-15.3

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.