

Infectious Disease Control IV– Level 4

Assignment of values(LOT 22K0701)

Instrument	Analyte	Unit	Infectious Disease Control IV-Level 4	
			Target	Range
Auto Lumo A2000/Auto Lumo A2000 plus	Anti-HBs	mIU/mL	213.27	149.29-277.24
	Anti-HBe	PEI U/mL	2.22	1.55-2.88
	Anti-HBc	PEI U/mL	20.27	14.19-26.36
	HIV p24 Ag	S/CO	12.18	8.52-15.83
Auto Lumo A1000	Anti-HBs	mIU/mL	220.67	154.47-286.87
	Anti-HBe	PEI U/mL	2.10	1.47-2.73
	Anti-HBc	PEI U/mL	20.75	14.53-26.98
	HIV p24 Ag	S/CO	11.94	8.36-15.52

Remark: The range printed in the assigned value table should only be considered as a reference and it is recommended that each laboratory establish its own mean and acceptable ranges.