Autobio Diagnostics Co., Ltd.

offers more than 600 in vitro diagnostic products including CLIA (microplate based CLIA and magnetic particle based CLIA), ELISA, POCT (Point of Care Test), Microbiology and Biochemistry. As an ISO9001 and EN ISO13485 manufacturer, Autobio supplies high quality products through its well established sales network and is renowned as a reliable partner. For details please visit en.autobio.com.cn. Autobio Diagnostics Co., Ltd. | NO.87 Jingbei Yi Rd | National Eco&Tech Zone | Zheng zhou City | China 450016.

You can find us on:
Facebook: @AutobioGlobal
Youtube Channel: Autobio Diagnostics
Email: info@autobio-diagnostics.com
WhatsApp: +8618595767835
Tel: +86 371 6200 7036







2025 AUTOBIO PRODUCT LIST





About Autobio

Autobio Diagnostics Co., Ltd was founded in 1998, headquartered in National Eco & Tech Zone, Zhengzhou.

Autobio specialized in R&D, production, marketing and service of clinical diagnostic products, especially covering immunoassay, microbiology, biochemistry coagulation and molecular fields, providing a comprehensive solution for medical laboratories. Autobio was listed on Shanghai Stock Exchange on September 1st, 2016, which was the first IVD manufacturer listed on the main board in Shanghai, China.

Autobio has become one of the leading companies with the most number of registered products and the most comprehensive product catalogue, currently having 794 products registration certificates, 612 CE certifications and 2342 patents.

Autobio has passed the ISO9001 and ISO13485 certifications, and strict quality management leads to continuous quality improvement.

Autobio's marketing network covers the whole country, and our products have entered more than ten thousands medical institutions. In oversea markets, Autobio's products has performed enormous potential, exporting to over 100 countries such as Europe, Asia, Middle East, South America, Africa, etc.

Contents

Autobio Product List


Immunosassay Solutions01
Reagents01
Instruments03
Clinical Chemistry Solutions05
Reagents05
Instruments
Microbial Solutions
Autobio Microbial Total Solution08
Reagents08
Instruments09
Molecular Solutions
Reagents12
Instruments13
Coagulation Solutions12
Reagents12
Instruments14

Immunoassay Solutions

Autobio

Reagents

CLIA Microparticles










療	Respiratory Disease			
M.	pneumoniae IgM	L. pneumophila IgM		
M.pneumoniae IgG Enterovirus 71 IgM				
C.	Coxsackie B IgM			
	RSV IgM	Influenza A IgM		
	Adenovirus IgM	Influenza B IgM		
	Adenovirus IgM	Influenza B IgM		

	Tuberculosis	
01	TB-IGRA	



Bone Meta	bolism
Osteocalcin Calcitonin	25-OH Vitamin D PTH



Hepato Fib	rosis	
PIIINP	НА	
Col IV	CG	
LN		

6	TORCH	
	Toxo IgM CMV IgM Rubella IgM HSV-1 IgM	Toxo IgG CMV IgG Rubella IgG HSV-1 IgG
	HSV-2 IgM	HSV-2 IgG

Cardiac Markers		
MYO	hs-cTnT	
cTnl	BNP	
CK-MB	H-FABP	
NT-proBNP		

Infectious Diseases				
Anti-HCV	HP-lgG			
HAV-IgM	Anti-HTLV*			
HEV-IgM	EBV VCA-IgA*			
HEV-IgG	EBV EA-IgA*			
Anti-HIV	EBV VCA-IgM*			
HIV Ag/Ab Combo	EBV NA1-lgA*			
Anti-TP	EBV EA-lgG*			

₩ Thyroid	
ТЗ	Anti-TPO
T4	TG
TSH	PTH
FT3	rT3
FT4	TRAb
Anti-TG	


6	Pre-natal	
	Free β-HCG P-AFP β-HCG	uE3 PAPP-A Inhibin A*

••	Anemia	
	Folate Vitamin B12	Ferritin
*	Rheumatoid Arth	ritis

RF IgG

	R	F IgM	
₩]	Allergy		
		tlgE	

Anti-CCP IgG





Antiphospholipid Syndrome

CLIA >>>



ELISA >>>



	Tumor Markers	TORCH	
	AFP HCG	Toxo-IgM	HSV-2 IgG
	CEA	CMV-lgM	Toxo-IgG
		Rubella-IgM	CMV-lgG
		HSV-1 IgM	Rubella-IgG
		HSV-2 IgM	HSV-1 IgG
*coming soon			

Endocrine Control II

Endocrine Control II Level 1
Endocrine Control II Level 2
Endocrine Control II Level 3

Cardiac Markers Control

Cardiac Markers Control Level 1
Cardiac Markers Control Level 2
Cardiac Markers Control Level 3

Tumor Marker Control II

Tumor Marker Control II Level 1
Tumor Marker Control II Level 2
Tumor Marker Control II Level 3
Tumor Marker Control II Level 4

ToRCH IgG Control

ToRCH IgG Control Level 1
ToRCH IgG Control Level 2
ToRCH IgG Control Level 3
ToRCH IgG Control Level 4

ToRCH IgM Control

ToRCH IgM Control Level 1
ToRCH IgM Control Level 2
ToRCH IgM Control Level 3
ToRCH IgM Control Level 4

Instruments

01

- · Small but fast
- Long walk away time
- · Low cost of consumbles

AutoLumo 5900



02

- Smarter: versatile and strong, brings more accurate results
- ·Simpler: easy operation and time-saving
- · Better: superior performance

AutoLumo A 1860



03

- Robust Performance
- · Easy Operation
- · Scalable Flexibility

AutoLumo #2000 Plus



04

- · Ultra-throughput
- · Flexible Scalability
- · Dual-core
- · Comprehensive Menu
- Convenient Experience

AutoLumo A6200



Clinical Chemistry Solutions

Autobio Product List

Instruments -



AutoChem 8801

• High Throughput

Constant 800 photometric tests/hour, 600 tests/hour for ISE

Advanced Optical System

Flat-field holographic concave grating photometer with 16 wavelengths

• Unified Operation Experience

Easy to operate with less training time


• Comprehensive Testing Panel



• AutoLumo A2000 Plus

AutoChem B2000 Series

Automatic Chemistry Analyzer



High-quality and accurate new experience, Biochemical testing more simple

Reagents*-

D Liver

5'-Nucleotidase

Adenosine Deaminase

Glutamate Dehydrogenase

Glycyl Prolyl Dipeptide Aminopeptidase

Glutathione Reductase

Leucine Aminopeptidase

Mitochondrial AST

Albumin

Alkaline Phosphatase

ALT (GPT)

AST (GOT)

Bilirubin Direct

Bilirubin Total

Ceruloplasmin

Cholinesterase

Gamma-GT

Prealbumin

Total Bile Acids

Total Protein

ULipids

Apo A1

Apo All

Аро В

Apo CII

Apo CIII

Apo E Cholesterol Total

HDL-C

LDL-C

Lipoprotein (a)

Phospholipid

Triglycerides

Electrolytes

Ammonia

Calcium

CO

Copper

Iron

Magnesium

Phosphorus

Zinc

🔀 Renal

Microalbumin

β-N-Acetylglucosaminidase

Neutrophil

Gelatinase-associated Lipocalin

Retinol-Binding Protein

Urine Transferrin

Cystatin C

Urine/CSF Total Protein

α 1-Microglobulin β 2-Microglobulin

Creatinine

Urea

Uric Acid

Diabetes

D-3-hydroxybutyrate Glycated Albumin

Non-Esterified Fatty Acids

Fructosamine Glucose

HbA1c

Lactate

Anemia

Unsaturated Iron Binding Capacity

G6PD Haptoglobin

Transferrin Ferritin

Pancreatitis

Amylase Pancreatic

Lipase

Other

Serum Amyloid A

Angiotensin Converting Enzyme Alpha-L-Fucosidase

Inflammatory

Cardiac

Lactate Dehydrogenase

Superoxide Dismutase

Lactate Dehydrogenase 1

Heart Fatty Acid Binding Protein

Lipoprotein-associated Phospholipase A2

Rheumatism & Immune

α-Hydroxybutyrate Dehydrogenase

Homocysteine

CK

CK-MB

hs-CRP

D-Dimer

Anti-Streptolysin O

Rheumatoid Factor

CRP

ΙgΑ

IgG

ΙgΜ

ΙgΕ

СЗ

C4


Myoglobin

*The products are still under registration and are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed.

Microbial Solutions

Autobio **Product List**

Autobio Microbial Information System



Reagents -

Mycoplasma ID/AST Reagents

Mycoplasma(UU-MH) Mycoplasma IES

Mycoplasma TIES

Mycoplasma IES Plus

Blood Culture System

Lytic Anaerobic Culture Bottle

Standard Aerobic/Anaerobic Culture Bottle

Resin Aerobic/Anaerobic Culture Bottle

Standard Pediatric Culture Bottle

Resin Pediatric Culture Bottle

Mycobacteria Culture Bottle

Q MALDI-TOF MS System

Sample Pretreatment Reagent

Calibrator for use with AUTOF MS

CHCA Matrix for use with AUTOF MS

Positive Blood Cultures Pretreatment Reagent

Filamentous Fungi Pretreatment Reagent

MALDI-TOF MS Control 5T/20T

ID/AST System

Gram-positive, GP Enterobacteriaceae, EB

Non-fermenters, NF Streptococcus, ST

Fungal, FG

Gram Negative Bacteria ID/AST Gram Positive Bacteria ID/AST Yeast ID/AST

Gram Negative Bacteria ID Gram Positive Bacteria ID

Instruments



AutoStreak \$1800

• High efficiency, treatment speed of 180

plates/hour

International patent inoculation brush,

10 μ L quantitative inoculation

Independent storage area of abnormal

samples

Real-time monitoring of consumption and

experiment progress

- Higher plate utilization

•••

Support the expansion of up to 360 bottle

positions

Higher flexibility for integrated control unit

Better temperature uniformity

• - Low noise (<45 dB) and small width



B(60



• 120 cell capacity

The instrument can carry out cascade connection

High-precision temperature control system keeps constant temperature of culture bottle

Continuous detection in each 10 minutes, minimum positive sample detection time is 3 hours



High throughout

0.1 second report time

Complete pretreatment reagents

Simple sample pretreatment


Faster detection speed

AST and Nucleic acid MS function



Autof **MS**





Autof

B(120 Plus 10

Molecular Solutions

Autobio

• Fully automated system, automatic sample addition, incubation, detection and result output.

High throughput, true MIC detection, real-time detection.

120-well card design, including more antimicrobial drugs and concentrations

Multiple card types: ID cards, AST cards, ID/AST combo cards.

Use the latest version of CLSI, EUCAST standards.

Provide complete fungal AST solutions.



AutoMic-1600

Pathogen in Respiratory Tract >>>

AutoMolec SARS-CoV-2 rRT-PcR

AutoMolec M. Pneumoniae/Adenovirus

AutoMolec M. Pneumoniae & Macrolides -Resistance

AutoMolec M. Pneumoniae/C. Pneumoniae/Adenovirus*

AutoMolec FluA/FluB/RSV

AutoMolec SA/MRSA/KPN

AutoMolec PAE/HIN/E.coli

AutoMolec Victoria/Yamagata

AutoMolec SP/LP/MC

AutoMolec AIV H5/H7

AutoMolec Flu A H1/H3

AutoMolec HPIV

AutoMolec HMPV

Reagents

AutoMolec HPV 16/18

Women's Health >>>

AutoMolec GBS

AutoMolec HSV 1/2

AutoMolec HPV 6/11

AutoMolec CT/NG/UU*

Individualized Therapy >>>

AutoMolec CYP2C19

AutoMolec ALDH2

AutoMolec MTHFR/MTRR

Gastrointestinal >>>

AutoMolec Norovirus

AutoMolec Rotavirus group A

AutoMolec Salmonella/Shigella

Infectious Diseases >>>

AutoMolec Monkeypox

AutoMolec HBV DNA Quant*

AutoMolec HCV RNA Quant*

AutoMolec HIV-1*

*RUO

Instruments



. . .

Small sample batch management

6 assays on-board at any one time

16 independent 4-well temperature control modules

Single sample multiple detection

Supporting both RNA and DNA

160 results in an 8 hour shift



AutoMolec 3000

Random access

Comprehensive anti-contamination solution

Multiple items are detected simultaneously

41 independent control channels

60 samples (circular rotor), priority lane for

emergency samples

20 items at one time

200 tests/8 hours

Coagulation Solutions

Autobio Product List

Reagents And Consumables -

- Prothrombin Time (PT) Reagent
- Activated Partial Thromboplastin Time (APTT)
- Thrombin Time (TT) Reagent
- Thrombin Reagent (for Fbg Assay)
- Imidazole Buffered Saline
- Coagulation Control Plasma
- D-Dimer Assay Kit

- D-Dimer Control
- Fibrin Degradation Product (FDP) Assay Kit
- Fibrin Degradation Product (FDP) Control
- Antithrombin III Assay Kit
- Antithrombin III Control
- Clean Solution
- Rinse Solution
- Reaction Cuvette

Instruments



Ci-310




• - Constant, 300 tests/hour

Multi-wavelength technology

Automated mixing study

Pack-and-place cuvette transport technology

Original liquid, ready-to-use reagents



Ci-120x

+- PT: 150 tests/hour

Multi-wavelength technology



Automated mixing study

Built-in sample and reagent barcode reader

Original liquid, ready-to-use reagents